

THE ONLY DEVICE

in the WORLD that FEELS your BODY

For non-surgical BODY SHAPING and SKIN PERFECTION



• MICROVIBRAZIONE COMPRESSIVA® ENCOSPHÈRES® THERAPY

ENDOSPHÈRES THERAPY® is a brand-new treatment harnessing the power of COMPRESSIVE MICROVIBRATION®, patented system engineered for non-invasive skin rejuvenation, body shaping and physical therapy. The method has been studied by influential medical research centers, that have certified the efficacy of the treatment.

The treatment is delivered using a cylinder made up of special 55 silicone spheres positioned in a honeycomb-like pattern, that have a specific density and diameter.


SENSOR


TECHNOLOGY

Our innovative SENSOR technology increases the efficacy of the Endosphères Therapy treatment.
How does this work?

The sensors positioned inside the delivery handpiece detect tissue resistance and they determine the correct pressure required to treat each area optimally. This means it is easier to carry out

INDIVIDUAL BESPOKE TREATMENTS and ensure best results.





"HONEYCOMB" arrangement



The device uses COMPRESSIVE MICROVIBRATION to transmit low frequency vibrations creating a pulsed rhythmic action on the tissue. The wave resulting from the vibrations has 5 therapeutic benefits on the tissue. It provides a well-rounded treatment. Initially it restores the normal state of tissues and microcirculation and in later stages works on reshaping and contouring the body.




VASCULAR action

LYMPHATIC drainage

body CONTOURING

MUSCLE toning



Compressive Microvibration acts on mechanoreceptors to reduce pain and tissue inflammation. It improves tissue oxygenation, reducing painful swelling.

The honeycomb arrangement of the spheres stimulates the vascular system, affecting the blood flow in the tissues. The "vascular workout" of the tissue improves microcirculation.

Compressive Microvibration uses the rotational movement of the spheres and the direction of the cylinder on the tissue to activate the lymphatic system to work correctly. This helps to eliminate excess fluids that have accumulated as a result of poor circulation.


The oscillations of the spheres create a push pull action which remodels the tissue. The shape and contour of the body is improved after only a few sessions, as well as reducing unsightly orange peel effect visible on skin.

Compressive Microvibration technology tones the muscles as they actively resist the compression. The vibration stimulates the muscles and can improve tone in the treated areas.






- PATENTED METHOD
 100% MADE IN ITALY
- AUTOMATIC BOOST WITH POWER UP TO 355Hz
- LAST GENERATION
 TECHNOLOGY
- PRESSURE SENSOR
- an **ADDITIONAL CYLINDER**may be requested to treat specific vascular conditions







15 shore gel



The low frequency mechanical movement of the spheres work deep into the tissue and act on the main causes of cellulite

BOOST 355 Hz

DYNAMIC frequency/pressure INDICATOR/

> REMOTE CONTROL JOYSTICK



Ak SENSOR all



FACIAL TREATMENT Areas of application

Mouth • Eyes • Forehead • Cheeks Cheekbone • Glabella Neck and lower neck area • Breast Pre and post Injectable treatments

- **SMOOTHES WRINKLES**
- MPROVES VASCULARIZATION AND OXIGENATION
- LIFTS AND FIRMS THE SKIN
- ▶ RELAXES CONTRACTED MUSCLES (EXPRESSION LINES)
- STIMULATES COLLAGEN PRODUCTION



COMFORTABLE for patients







BODY TREATMENT Areas of application

Cellulite • Orange Peel skin • Contracted Muscles • Oedema Swollen legs • Arms • Buttocks • Ankles • Abdominal fat Breast tightening • Post-pregnancy treatment

- REDUCES LYMPHATIC FLUID RETENTION
- REDUCES ORANGE PEEL SKIN
- MPROVES VENOUS CIRCULATION
- IMPROVES BODY SHAPE
- TONES MUSCLES



INCREASE OF salon/clinic revenue



The FIRST device in the world to AUTOMATICALLY change the intensity of the microvibrations according to tissue resistance.



Clinical observations have confirmed Endosphères Therapy as a valid and trustworthy treatment providing a solution for cellulite, poor venous circulation, fluid retention and muscle tension.

In the studies, in order to obtain the correct diagnosis, patients were subjected to a clinical evaluation. This included the inspection and palpation of the tissue. Echo colour Doppler, echography, thermography and biometrical impedance analysis were carried out. Analysis was carried out after every six sessions to give an insight into the tissue changes. The results show the efficacy of Compressive Microvibration.

CLINICALLY PROVEN RESULTS











Results after 18 sessions













Results after 3 sessions





Results after 3 sessions

Results after 3 sessions















Ak Sensor All

TECHNICAL DATA

Diameter

Height

Weight

Vibration

Cooling

Sensor System

Touch Screen

Sound Level

Electrical Features

315mm(up)/364mm(down)

1.150 mm

18 kg

29 - 355 Hz

ED/AIRCOOLING

FND-21K

Capacitivo 7 Colour

56dBA

da 100 a 240 V./50-60Hz

www.endospheres.com

The Endosphères Therapy Ak Sensor All patent belongs to the Fenix Group and is manufactured by the Fenix Group



The Endosphères Therapy brand belongs to the Fenix Group, a 100% Italian company. The group has been at the forefront in the field of aesthetics and medical-aesthetic equipment, always basing its strength on Research and development.

The Fenix Group has developed Endosphères Therapy, a new technology, using a patented Compressive Microvibration method. In recent years this has proved itself as the treatment of excellence to combat edematous fibrosclerotic panniculopathy (cellulite), adipose accumulations and face rejuvenation.

The Fenix Group is an innovative company investing in

marketing and PR campaigns both in trade press and consumer press. It is active across all social media platforms. Customers are assisted with in salon/clinic marketing to promote the treatments locally in their area. Service and technical support is available throughout the world.

The Fenix Group's mission is to design high quality electro-medical devices and manufacture the most innovative, effective and science-based technologies. These essential features ensure safety and efficacy making the Fenix Group a successful and customer-driven company.



Headquarters ITALY

Fenix Group s.r.l. - Strada Lungofino, 187 - c/o C.C. Ibisco B9/10 Città Sant'Angelo (PE) Tel. +39 085 4483235 - info@fenixgroup.it



Branch Office U.K.

Fenix Group LTD - 51 Holland Street, London W8 7JB Tel: +44(0)2070960040 - ukmanagement@fenixgroup.it



Branch Office SPAIN

Fenix Endospheres Therapy SL - C/ Cirilo Amoros, 6 Entresuelo 46004 Valencia Tel. +34 627 508380 - espmanagement@fenixgroup.it



Branch Office FRANCE

Fenix France S.A.S. - Startway 43-47 Avenue de la Grande Armee - 75016 PARIS Tel: +33(0)170395146 - frmanagement@fenixgroup.it

For other countries look for distributors on www.endospheres.com